# AXCENT ACRYLIC™

## Orthodontic Resin

### CADMIUM FREE

### Processing Instructions

Complies with ISO 20795-2:2010 Specification Type I, Class 1

Use Axcent Acrylics Orthodontic White and Black Powders with Axcent Acrylics #3 Self-Cure Liquid. Use Axcent Acrylics Orthodontic Clear with Self-Cure Pink and #4 Hi-Clarity Monomer Liquid.

#### INTENDED USE:

Fabrication of orthodontic appliances.

#### CAUTION:

This product is to be used only under the supervision of a licensed dental practitioner.

#### INDICATION(S) FOR USE:

Garreco Axcent Acrylics Orthodontic Resin is a self-cure acrylic resin used for orthodontic applications.

#### WARNINGS AND PRECAUTIONS:

Liquid contains methyl methacrylate monomer. Highly flammable, keep away from sources of ignition. NO

SMOKING. Irritating to eyes, skin and respiratory system. IF ON SKIN: May be irritating to skin in some sensitive individuals, especially after prolonged and/or repeated contact. IF INHALED: May be irritated by gross overexposure, no matter how generated. IF SWALLOWED: Not expected to cause any harmful effects. Ingestion of large amounts may be irritating if product is swallowed, may cause nausea, headache, vomiting and/or diarrhea. IF IN EYES: May be irritating to the eyes by gross overexposure, no matter how generated. Symptoms of overexposure may include redness, itching, irritation and watering. Keep dust out of eyes.

#### CONTRAINDICATION(S):

Do not use with patients with a known allergy to methacrylates or other components of the product.

#### ADVERSE REACTION(S):

Reversible inflammatory changes of the oral mucosa may occur with susceptible persons.

#### INSTRUCTIONS:

If model is dry, submerge in water until bubbling ceases; this is to ensure all air is removed from the stone. Coat surface with isolating solution and allow to dry thoroughly. Assemble wire framework on model.

Application -

The "salt & pepper" technique is recommended. First, wet surface of model with liquid, then sprinkle powder on model in a thin layer. Saturate powder by slowly adding liquid with a dropper bottle. Continue adding powder and liquid alternately until the desired thickness is achieved. Too much liquid added at one time may result in excessive flow, which can be controlled by the addition of more powder. Trim excess material while still in rubbery dough stage by using a sharp instrument. Smooth surface of appliance using liquid and a brush.

Allow time to cure or use a hydroflask to prevent formation of bubbles. Cure time is approximately 20 minutes at 23°C (73°F) or 10 minutes in hydro flask at 40-45°C (104-113°F) and 35-40 lbs. pressure.

Trim and polish in normal manner.

#### STORAGE:

Dental acrylics based on methyl methacrylate should be stored in a cool (16-26°C, 61-79°F), dry place. Avoid prolonged exposure to sunlight. Keep containers closed when not in use.

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Garreco, LLC