AXCENT ACRYLIC™

Dynasty Hi-Impact Acrylic Denture Base

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CADMIUM FREE

Processing Instructions

Complies with ISO 20795-1:2008 Specification Type I, Class 1

Use with Garreco Axcnet Acrylics Heat Cure Acrylic Liquid Monomer

INTENDED USE:

For the fabrication of dental appliances.

CAUTION:

This product is to be used only under the supervision of a licensed dental practitioner.

INDICATION(S) FOR USE:

Garreco Axcnet Acrylics Dynasty Acrylic Denture Base are Heat-cure and microwaveable acrylic resin for denture bases with a long working time, strength and flexibility.

WARNINGS AND PRECAUTIONS :

LIQUID MONOMER - Highly flammable liquid and vapor. IF INHALED: May cause respiratory irritation. IF ON SKIN: May cause skin irritation and allergic skin reaction.

POWDER - IF ON SKIN: May be irritating to skin in some sensitive individuals, especially after prolonged and/or repeated contact. IF INHALED: May be irritated by gross overexposure, no matter how generated. IF SWALLOWED: Not expected to cause any harmful effects. Ingestion of large amounts may be irritating if product is swallowed, may cause nausea, headache, vomiting and/or diarrhea. IF IN EYES: May be irritating to the eyes by gross overexposure, no matter how generated. Symptoms of overexposure may include redness, itching, irritation and watering. Keep dust out of eyes.

CONTRAINDICATION(S):

Do not use with patients with a known allergy to methacrylates or other components of the product.

ADVERSE REACTION(S):

Reversible inflammatory changes of the oral mucosa may occur with susceptible persons.

RECOMMENDED STORAGE CONDITIONS:

Store in original container. Store containers in a cool, dry location, away from direct sunlight, other light sources, or sources of intense heat. Keep container closed to prevent water absorption and contamination ignition and light.

DISPOSAL:

LIQUID MONOMER - This has been classified as a hazardous waste by the EPA under RCRA. After addition of excess inhibitor, dispose waste material in accordance with Federal, State and Local regulations.

POWDER - Dispose in a landfill or incinerate according to Federal, State, and Local regulations.

STANDARD INSTRUCTIONS:

No Trial Pack Method

1. Wax, invest, and boil out the case as usual.
2. Apply Separating Film to the gypsum surfaces before mixing the acrylic.

DFU052 rev.170322

3. Mixing and packing:

- a. Add by volume 3 parts powder to 1 part of #1 heat cure liquid (23.4 g powder / 10 ml liquid), spatulate approximately 45 seconds. Let stand until mix is of doughy consistency (5-8 minutes). Transfer to a moist warm flask, not more than 45°C (115°F) and pat into place with fingers covered with moist cellophane or latex gloves. Be certain the flask contains an excess of material before closing.
- b. Close flask immediately to full metal closure. Overflow should be observed on complete perimeter of flask. If overflow is noted, the flask is full and need not be reopened.
- c. Let flask bench set for 15 to 30 minutes before beginning the curing process. This allows the acrylic to reach a complete equilibrium.
- d. Cure at 73°C±1°C (163.4°F±1.8°F) for 8 hours; or 73°C±1°C (163.4°F±1.8°F) for 1 1/2 hours and then for 30 minutes in boiling water. Bench cool approximately 45 minutes before deflasking.

Standard Trial Pack Method

Use the same proportions of powder and liquid (3 to 1 by volume). Do not pack until the mix has a doughy consistency (5-8 minutes at 23°C/±1°C (73.4°F±1.8°F). Trial pack and cure as in Step 3d above.

Repairs

Repairs to processed Standard Denture Base material can be made with Denture Base Repair Material, following the instructions provided. Alternatively, the repair may be processed using the original Standard Denture Base material from which the denture was processed. Pack and cure the repair according to Standard Denture Base techniques.

MICROWAVE INSTRUCTIONS

CAUTION: Only use a microwave safe flask with this method.

1. Consult microwave owner's manual to determine the correct settings to achieve the following wattage levels:

Low setting 86-90 watts. High setting not more than 500 watts.

2. Packing

- a. Mix acrylic powder, 3 parts powder (by volume) to 1 part liquid (by volume).
- b. Allow to stand until resin has reached doughy consistency.
- c. Pack and/or trial pack as usual.
- d. Allow to bench set (compressed) for 30 minutes.

3. Use the appropriate steps below according to microwave configuration

Option 1 - Curing (Option 1 with Turntable and Wattage Control)

- a. Place flask vertically on turntable center. Cure by microwaving on low setting for 13 minutes.
- b. Reset microwave to high setting for 1 minute.
- c. Bench cool for 30 minutes, deflask and finish as usual.

Option 2 - Curing (Option 2 without Turntable but with Wattage Control)

- a. Place flask flat in the middle of microwave. Cure by microwaving on low setting for 6 ½ minutes. b. Turn flask over and cure an additional 6 ½ minutes (low setting).
- c. Turn flask over and change to high setting, cure for 30 seconds.
- d. Turn flask over and cure an additional 30 seconds (high setting).
- e. Bench cool for 30 minutes, deflask and finish as usual.

Option 3 - Curing (with Turntable but without Wattage Control)

- a. Place a cup of water inside with flask vertical and then cure by microwaving for 2 minutes.
- b. Rest for 5 minutes (in microwave).
- c. Cure for 5 minutes.
- d. Bench cool for 30 minutes, deflask and finish as usual.

Option 4 - Curing (without Turntable and without Wattage Control):

- a. Place a cup of water inside with flask flat and then cure by microwaving for 1 minute.
- b. Turn flask over and cure for 1 additional minute.
- c. Rest for 5 minutes (in microwave).
- d. Cure for 2 1/2 minutes.
- e. Turn flask over and cure for 2 1/2 additional minutes.
- f. Bench cool for 30 minutes, deflask and finish as usual.

Additional Notes:

If "bleaching" should occur, reduce high setting (wattage).

Appliances containing a greater or lesser acrylic volume may require adjustments to curing times.

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Garreco, LLC