



Garreco

SELF CURE

AXCENT ACRYLIC TRAY MATERIAL FOR CUSTOM TRAYS

Tray Material by Garreco is auto curing and designed for custom trays.

- Self curing acrylic
- Blue or white color available

Indication for use: Custom dental trays

CURE TYPE Self Cure
 INJECTABLE No
 MICROWAVABLE No
 POWDER/LIQUID RATIO 2:1



CUSTOM TRAY MATERIAL COLORS AND SIZES AVAILABLE

BLUE	WHITE
1 LB JAR - 3301501S	1 LB JAR - 3303001S
1 LB CASE (4x 1LB) - 3301501	1 LB CASE (4x1 LB) - 3303001
5 LB JAR - 3301505S	5 LB JAR - 3303005S
5 LB CASE (4x 5LB) - 3301505	5 LB CASE (4x5 LB) - 3303005
20 LB PAIL - 3301520	20 LB PAIL - 3303020



MATCHING LIQUID MONOMER

Liquid Monomer #6 - Tray Liquid
 8 OZ BOTTLE - 3920008
 8 OZ CASE (4x8 OZ) - 3920008S
 1 QT BOTTLE - 3920030
 1 QT CASE (4x1 QT) - 3920030S

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com