



Garreco

DIE STONE

CAD-SCAN TYPE IV-R RESIN DIE STONE

CAD-scan is a unique gypsum material that contains a special scanning additive and resin reinforcement. CAD-scan reflects light nicely for essential readability, allowing for a shorter working time, and is manufactured from an ultra-fine gypsum base, which prevents nicks and abrasions caused by brittleness. CAD-scan has a low expansion and high compressive strength, making it the material of choice for all CAD/CAM restorations. It is perfect for production environments, as well as those single cases requiring a rapid set.

- Excellent surface scanning
- Ideal for CAD/CAM restorations
- Resin reinforced
- Low expansion (0.06-0.09%) for accuracy

Indication for use: CAD/CAM Restorations

SETTING TIME 11-13 min
 WATER/POWDER RATIO 20ml:100g
 COMPRESSIVE STRENGTH 18,500 psi
 EXPANSION 0.06-0.09%



GOLDEN GREEN BUFF



10LB JAR	1466710	1462510	1461010
25LB PAIL	1466725	1462525	1461025
50LB	1466750	1462550	1461050

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com