



# Garreco

## PLASTER

### DENTAL PLASTER TYPE II PLASTER

Our Dental Plaster is formulated especially for use in today's complex field of dentistry. Designed for articulating and flasking procedures, Garreco dental plaster mixes smoothly and exhibits excellent working characteristics. Available in three formulations with varying set times.

- Slow, regular, or fast set times
- Smooth mixing
- Low expansion (0.19%)
- Excellent working characteristics

**Indication for use:** Articulating and flasking procedures

SETTING TIME ..... 17-19 minutes (Slow);  
 13-15 minutes (Medium); 9-11 minutes (Fast)  
 WATER/POWDER RATIO ..... 45ml:100g  
 COMPRESSIVE STRENGTH ..... 3,500 psi  
 EXPANSION ..... 0.19%



WHITE



25LB SLOW	1550025
50LB SLOW	1550050
25LB REG	1450025
50LB REG	1450050
25LB FAST	1500025
50LB FAST	1500050

Manufactured by Garreco, LLC  
 430 Hiram Road, Heber Springs, AR 72543, USA  
 Toll Free: (800) 334-1443  
 garreco.com