



Garreco

PERFECT PRESS INVESTMENT MATERIAL

Perfect Press is a graphite-free, phosphate bonded precision investment developed especially for pressed ceramics. Utilizing a fast casting technique, Perfect Press not only divests easily, but delivers accurate consistent fits. Conventional or fast casting. For use with Hi-X or Regular investment/expansion liquids.

- All ceramic systems
- All pressable systems
- Carbon free
- Convenient pouches
- Divests easily

Indication for use: Investment for pressed ceramics

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com



PERFECT PRESS SIZES AVAILABLE

- 50x100 G POUCHES - 6150016
- 100x100 G POUCHES - 6150026




MATCHING INVESTMENT LIQUID

- 1 L BOTTLE - 6200003S
- 1 L CASE (4x1 L) - 6200412
- 1 L CASE (16x1 L) - 6201647

Requires 1 liter of expansion liquid per 5 kilograms of investment material

LIQUID REQUIREMENTS

- 50x100 G - 1 LITER
- 100x100 G - 2 LITERS

*  **WARNING: Cancer** - www.P65Warnings.ca.gov

Garreco, LLC is an ISO 13485 Facility
Garreco is a trademark of Garreco, LLC. Heber Springs, AR. ©2019 All rights reserved.