



Garreco

SPECIALTY STONE

X-ROC EXPANSION STONE

X-Roc was developed for pouring the duplicate model in processing thermoplastic appliances such as: Flexite®, Valplast®, Flex-Star V, TCS®, FRS®, DuraFlex®, and DurAcetal®. This stone compensates for the shrinkage of these materials, thereby eliminating tight fits on removable prosthetics processed with conventional stones.

- Eliminates tight fits
- Moderate PSI
- Easy to work with

Indication for use: Specialty duplication material

SETTING TIME 22-25 min
 WATER/POWDER RATIO 29ml:100g
 COMPRESSIVE STRENGTH 13,000 psi
 EXPANSION 1.10%



WHITE



150G CASE (15x150G)	1483005
25LB	1483025
50LB	1483050

Manufactured by Garreco, LLC
 430 Hiram Road, Heber Springs, AR 72543, USA
 Toll Free: (800) 334-1443
 garreco.com