




# Garreco



## BASE PLATE WAX



### MEDIUM SOFT

#### Base Plate Wax

Garreco Medium Soft Base Plate Wax is a dependable wax that boils out clean, is competitively priced, and polishes to a high shine.

**Sizes:** 76mm x 152mm x .060 thick, available in 1lb or 5lb

MEDIUM SOFT 1LB - 51510101

MEDIUM SOFT 5LB - 51510105

### MEDIUM HARD

#### Base Plate Wax

Garreco Medium Hard Base Plate Wax has the same characteristics as the Medium Soft, but is formulated with higher melting waxes and is specifically designed for warmer climates.

**Sizes:** 76mm x 152mm x .060 thick, available in 1lb or 5lb

MEDIUM HARD 1LB - 51510201

MEDIUM HARD 5LB - 51510205

### ALL SEASON

#### Base Plate Wax

Our All Season Base Plate Wax was designed to be used most anywhere a good quality, economy wax is needed. This wax boils out clean, polishes nicely, and is formulated to securely hold teeth.

**Sizes:** 76mm x 152mm x .060 thick, available in 1lb or 5lb

ALL SEASON 1LB - 51510001

ALL SEASON 5LB - 51510005