



GINGIVA 3D - COMPATIBLE PRINTERS

Printer Name	Compatible?
MakeX M-One Pro 70	Yes
Miicraft Prime	Yes
Miicraft Ultra	Yes
Asiga Max UV	Yes
Asiga Pro 4K	Yes
Ackuretta Dentiq	Yes
Ackuretta Freeshape 120	Yes
Shining3D AccuFab-L4D	Yes
Phrozen Sonic Mini 4K	Yes
Phrozen Series	Preliminary settings
Peopoly Phenom Prime	Yes
Anycubic Photon	Yes
Anycubic Photon Mono	Yes
Anycubic Photon Mono SE	Yes
Anycubic Photon Mono X	Yes
Crealty Halot-One	Yes

Printer Name	Compatible?
Elegoo Mars	Yes
Elegoo Mars 2	Yes
Elegoo Mars 3	Yes
Elegoo Saturn	Yes
Moonray S & D	Not compatible
Formlabs Form2	Not compatible
385 – 420nm MSLA, SLA & DLP printers	Contact us