



## MODEL 3D ULTRA

Model 3D Ultra from Garreco Print is perfect for manufacturing orthodontic models and aligner models. This skin/peach colored resin is suitable for open high speed printing on Liquid Crystal Display (LCD) and Digital Light Processing (DLP) 3D-printers.

### PRODUCT DESCRIPTION

Model 3D Ultra is an opaque skin/peach colored photopolymer. Its color accentuates depth in your dental models, which accommodates a perfect view on undercuts in 3D-printed models. 3D-printed parts from this material have great dimensional stability and low shrinkage during printing. Model 3D Ultra is easy to use on all open source LCD and DLP 3D-printers in the range of 385 –420nm. This material has excellent properties like low shrinkage, low odor and high print speed making it perfect for the production of orthodontic and aligner models.

### KEY BENEFITS

- Good dimensional stability
- Great detail
- Temperature resistant for aligner production
- Low shrinkage

### 3D-PRINTER COMPATIBILITY

- Asiga Max and 4K
- Shining3d AccuFab-D1
- Photon and Phrozen Series
- All low power, open source 385 - 420nm LCD and DLP 3D-printers

### PRODUCT CONFIGURATIONS

Model 3D Ultra - 1 kg (1000 g)  
Item 6865502S



Distributed by Garreco, LLC  
430 Hiram Road, Heber Springs, AR 72543, USA  
Toll Free: (800) 334-1443  
garreco.com

Garreco, LLC is an ISO 13485 Facility  
Garreco is a trademark of Garreco, LLC. Heber Springs, AR. ©2023 All rights reserved.

Rev. 230213