



ULTRA-SONIC CLEANER #5 PLASTER AND STONE REMOVER

Ultra-Sonic Cleaner #5 is especially formulated to remove plaster and stone from dentures, bridges, crowns, orthodontic appliances, and other dental restorations.

- Removes plaster and stone
- Fast working

Indication for use: Removes plaster and stone from dental appliances

Standard Size Available: 1 Gallon (3.8 Liters)

Directions: Use full strength; do not dilute. Keep covered when not in use. **DO NOT USE ON ALUMINUM.**

Item #: 3610040



RELATED PRODUCTS

- Ultra-Sonic Cleaner #1: General Purpose Concentrate - 3570040
- Ultra-Sonic Cleaner #2: Tartar and Stain Remover - 3590040

Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com