



# Garreco

## NATURAL ROCK

TYPE IV-R RESIN REINFORCED DIE STONE

Natural Rock is a traditional, creamy, smooth, resin reinforced die stone which exhibits high compressive strength and very low expansion.

—Exceptionally dense for maximum strength

—Very low expansion (0.08%)

**Indication for use:** Crown and bridge; porcelain; implants; cast partials; high-precision cases

SETTING TIME .....	10-13 min
WATER/POWDER RATIO .....	20ml:100g
COMPRESSIVE STRENGTH .....	18,000 psi
EXPANSION .....	0.08%

LT. TAN



25LB	1446725
47.5LB	1446750

Manufactured by Garreco, LLC  
430 Hiram Road, Heber Springs, AR 72543, USA  
Toll Free: (800) 334-1443  
garreco.com

## DIE STONE



Check out our economically-priced alternative to FujiRock® Die Stone!