



PLASTER

ARTICULATING PLASTER

TYPE II PLASTER

Articulating Plaster exhibits similar characteristics to our Dental Plasters, but with a greatly reduced set time. Recommended for articulation use.

- Smooth mixing
- Low expansion (0.19%)
- Excellent working characteristics

Indication for use: Articulation

SETTING TIME 3-5 min.
WATER/POWDER RATIO 45ml:100g
COMPRESSIVE STRENGTH 3,500 psi
EXPANSION 0.19%

WHITE



25LB 1560025
50LB 1560050



Manufactured by Garreco, LLC
430 Hiram Road, Heber Springs, AR 72543, USA
Toll Free: (800) 334-1443
garreco.com