

Safety Data Sheet

Form No. A896

Date Prepared:

5/30/2023

SECTION 1 : Identification of the substance/preparation and of the company / undertaking**(a) GHS product identifier**

Garreco Silaform (All Types and Parts)

(e) Emergency phone numberCHEMTRAC 1-800-424-9300
CCN9105**(b) Other means of identification**

NA

(c) Recommended use of the chemical and restrictions on use

For professional dental applications.

(d) Supplier's detailsGarreco, LLC.
430 Hiram Road
Heber Springs, AR 72543
Phone: 1-800-334-1443**SECTION 2: Hazards identification****(a) GHS classification of the substance/mixture****Substance Name**

1. Polyorganosiloxanes with inert mineral fillers.
2. Platinum (Below Threshold)

(b) Label Elements**Hazard statements**

None

Precautionary statements**Hazard Symbol(s)**

NONE

Signal Word(s)

NONE

(c) Other hazards which do not result in classification

Slightly irritating to eyes

Combustible liquid. Explosive vapor/air mixtures may be formed above 200°C.

SECTION 3: Composition/information on ingredients**(a) Chemical(s) Identity:****(b) Common Name:**Polyorganosiloxanes with inert mineral fillers.
Platinum**(c) CAS No.**CBI
7440-06-4**Mixture:****Concentration (Percentage)**95.0 - 100.0%
0.004% (Below Threshold)**SECTION 4: First-aid measures****(a) Description of first aid measures:****IF ON SKIN (or hair):** Wipe off with cloth or paper. Then wash with soap and water.**IF INHALED:** No special measures required.**IF SWALLOWED:** In cases of sickness seek medical advice (show label if possible).**IF IN EYES:** Rinse immediately with plenty of water. Seek medical advice in case of continuous irritation**(b) Most important symptoms and effects, both acute and delayed:**

ND

(c) Indication of any immediate medical attention and special treatment needed:**IN EYES:** Tears, blurred vision and redness.**INGESTION:** Burning sensation of the mouth, throat and gastrointestinal tract and abdominal pain.**SECTION 5: Fire-fighting measures****(a) Suitable extinguishing media:**Water spray, foam, powders, carbon dioxide (CO²)**(b) Special hazards arising from the chemical or mixture:**

Combustible liquid. On heating or during combustion: Toxic and flammable vapors are released (explosion hazard).

(c) Special protective equipment and precautions for fire-fighters:

Self-contained breathing apparatus.

SECTION 6: Accidental release measures**(a) Personal precautions, protective equipment and emergency procedures:**

Safety Glasses.

(b) Environmental precautions:

This product does not present any particular risk for the environment.

(c) Methods and material for containment and cleaning up:

Absorb mechanically and dispose of according to regulations. Absorb with a liquid binding material such as diatomaceous earth and dispose of according to regulations.

SECTION 7: Handling and storage**(a) Precautions for safe handling:**

Use nitrile gloves only. Safety glasses.

(b) Conditions for safe storage, including any incompatibilities:

Stable under normal storage conditions. Incompatible with strong oxidizing agents.

SECTION 8: Exposure controls/Personal protection**(a) Control parameters:**

Chemical	ACGIH	OSHA
	TLV	PEL TWA
Platinum	1 mg/m	ND

(b) Appropriate Engineering Controls:

NA

(c) Individual protection measures:

Use safety glasses. Do not drink, eat or smoke in the workplace.

SECTION 9: Physical and chemical properties

(a) Appearance:	Viscous paste. Part A is grey/white. Part B is blue.
(b) Odor:	Odorless
(c) Odor threshold:	ND
(d) pH:	ND
(e) Melting point / freezing point:	ND
(f) Initial boiling point and boiling range:	ND
(g) Flash point:	>200° C
(h) Evaporation rate (BuAc=1):	ND
(i) Flammability:	ND
(j) Upper/lower flammability or explosive limits:	ND
(k) Vapor Pressure:	<0.01 kPa at 20° C
(l) Vapor density:	ND
(m) Relative density:	ND
(n) Solubility:	Practically insoluble in water. Dispersible in: diethylether, chlorinated solvents, aromatic or aliphatic hydrocarbons
(o) Partition coefficient: n-octanol/water:	ND
(p) Auto-ignition temperature:	ND
(q) Decomposition temperature:	>200° C
(r) Viscosity:	ND

SECTION 10: Stability and reactivity

(a) Reactivity:	Possible
(b) Chemical stability:	Stable at room temperature.
(c) Possibility of hazardous reactions:	Reacts with: chemical compounds with mobile hydrogen, metal salts and complexes. A fire hazard arises since highly flammable gas is released when in contact with alkalis, caustic products, strong oxidizing agents.
(d) Conditions to avoid:	Alkalis, caustic products, strong oxidizing agents.
(f) Hazardous decomposition products:	On combustion or on thermal decomposition releases: flammable vapors. Toxic gases (Carbon oxides). Quantity of hydrogen potentially released: 5 l/kg of product

Safety Data Sheet

Form No. A896

Date Prepared:**5/30/2023****SECTION 11: Toxicological information**

Acute toxicity	NA
Skin corrosion/irritation	ND
Serious Eye Damage / Irritation	May cause temporary irritation to ocular mucous membranes.
Respiratory or skin sensitization	ND
Germ cell mutagenicity	ND
Carcinogenicity	ND
Reproductive toxicity	ND
STOT-single exposure	ND
STOT-repeated exposure	ND
Aspiration Hazard	ND

(a) Exposure route: Eyes

(b) Symptoms related to the physical, chemical and toxicological characteristics:

May cause temporary irritation to ocular mucous membranes.

(c) Delayed and immediate effects and also chronic effects from short and long term exposure:

May cause temporary irritation to ocular mucous membranes.

(d) Numerical measures of toxicity:

LD 50 skin (Rat): >2000mg/kg.

LD 50 oral (Rat): >2000mg/kg.

SECTION 12: Ecological information

(a) Ecotoxicity:

Low to none

(b) Persistence and degradability:

Not Biodegradable

(c) Bioaccumulative potential

Not bioaccumulable

(d) Mobility in soil:

The product is slightly soluble and forms sediment by absorption in materials/particles in suspension.

(e) Other adverse effects:

ND

SECTION 13: Disposal considerations**Product:****Recommendation**

Do not discharge waste into drains. Dilute with a flammable solvent and incinerate at a licensed installation.

Un-vulcanized product and Vulcanized product: May be disposed of with non hazardous industrial waste.

Contaminated packaging: Allow it to drain thoroughly. Rinse with an appropriate solvent. Recover the solvent used for rinsing and incinerate at a licensed site. Entrust the decontaminated packaging to a licensed waste-contractor. Follow local, state, and federal regulations.

SECTION 14: Transport information

(a) UN Number

Not regulated.

(b) UN Proper shipping name

NA

(c) Transport hazard class(es)

NA

(d) Packing Group

NA

(e) Environmental hazards

NA

(f) Transport in bulk

NA

(g) Other Information

NA

Safety Data Sheet

Form No. A896

Date Prepared:**5/30/2023****SECTION 15: Regulatory information****SARA Reporting Requirements:**

Not required

SARA Threshold Planning Quantity:

NA

TSCA Inventory Status:

NA

Other Federal Requirements:

NA

Other Canadian Regulations:

NA

State Regulatory Information:

NA

SECTION 16: Other information**PREPARED BY:** Kristofer Mainar**GAR QMS SDS REFERENCE:** A040 and A041**HAZARDOUS MATERIAL IDENTIFICATION (HMIS) RATING:**

Health	ND
Flammability	ND
Reactivity	ND
Other	ND

REVISION NUMBER: 230530**CHANGES FROM PREVIOUS VERSION:**

Reviewed for accuracy.

ABBREVIATIONS

NA Not Applicable

ND Not Determined

NE Not Established

ppm parts per million

G Gallon

mg Milligram

L Liter

gm Gram

mol Mole

kg Kilogram

μ Micro

mm Millimeter

p Pico

Pa Pascals

c cent

LC Lethal Concentration

ACGIH American Conference of Governmental Industrial Hygienist

CPR Controlled Product's Regulation

OSHA Occupational Safety and Health Administration

PEL Permissible Exposure Limit

TLV Threshold Limit Value

LD Lethal Dose

TC Toxic Concentration

TD Toxic Dose

BOD Biological Oxygen Demand

COD Chemical Oxygen Demand

Lo Lowest

ThOD Theoretical Oxygen Demand

TLm Threshold Limit

IC Inhibitory Concentration

DOC Dissolved Organic Carbon

H Hours

M Months

D Days

Y Years

W Weeks

DSL Canadian Domestic Substances List

NDSL Canadian Non-domestic Substance List

IARC International Agency for Research for Cancer

NOEL No Observed Effect Level

NOAEL No Observed Adverse Effect Level

CBI Confidential Business Information

Safety Data Sheet

Form No. A896

Date Prepared: 5/30/2023

THIS MATERIAL SAFETY DATA SHEET IS PREPARED IN COMPLIANCE WITH FEDERAL REGULATIONS (29 CFR 1910.1200) OF CHEMICALS AND THE GLOBALLY HARMONIZED SYSTEM OF CLASSIFICATION AND LABELLING REVISION 5. ANY APPLICABLE STATE AND LOCAL REGULATIONS SHOULD BE CONSULTED.

THE ABOVE INFORMATION MAY BE BASED IN PART ON INFORMATION PROVIDED BY COMPONENT SUPPLIERS AND IS BELIEVED TO BE CORRECT AS OF THE DATE HEREOF. HOWEVER, NO WARRANTY OF MERCHANTABILITY, FITNESS FOR ANY USE, OR ANY OTHER WARRANTY IS EXPRESSED OR IS TO BE IMPLIED REGARDING THE ACCURACY OF THESE DATA, THE RESULTS TO BE OBTAINED FROM THE USE OF THE MATERIAL, OR THE HAZARDS CONNECTED WITH SUCH USE. SINCE THE INFORMATION CONTAINED HEREIN MAY BE APPLIED UNDER CONDITIONS BEYOND OUR CONTROL AND WITH WHICH WE MAY BE UNFAMILIAR, AND SINCE DATA MADE AVAILABLE SUBSEQUENT TO THE DATE HEREOF MAY SUGGEST MODIFICATION OF THE INFORMATION, WE ASSUME NO RESPONSIBILITY FOR THE RESULT OF ITS USE. THIS INFORMATION AND MATERIAL IS FURNISHED ON THE CONDITION THAT THE PERSON RECEIVING IT SHALL MAKE HIS/HER OWN DETERMINATION AS TO THE SUITABILITY OF THE MATERIAL FOR HIS/HER PARTICULAR PURPOSE AND ON THE CONDITION THAT HE/SHE ASSUME THE RISK OF HIS/HER USE THEREOF.